

Select the instruments you need to perform a complete diagnostic examination:



**Green Series™ 777 Wall Transformer**  
Energy efficient with handle-based rheostats and automatic on/off sensors make operation simple.



**3.5 V Coaxial Ophthalmoscope**  
Multiple aperture/filter options provide long-lasting, reliable performance for general and specialist examinations.



**PanOptic™ Ophthalmoscope**  
5x larger field of view of the fundus vs. standard ophthalmoscopes in an undilated eye (25° vs. 5°).



**MacroView™ Otoscope**  
A nearly complete view of the tympanic membrane with 30% greater magnification than a traditional otoscope.



**Pneumatic Otoscope**  
Rotating lens with ultraseal enables pneumatic otoscopy and easier instrumentation.



**KleenSpec® Disposable Specula Dispenser**  
Offers optimum organization of disposable otoscope tips. Comes filled with specula.



**Wall Aneroid**  
Easy-to-read, jeweled-movement, mercury-free manometer is certified to be accurate within one percent of scale ( $\pm 3$  mmHg).



**Spot Vital Signs®/Spot Vital Signs® LXi/  
Connex® ProBP™ Devices**  
Fast, accurate, and efficient way to take vitals. Allows for electronic transfer of data to EHR.



**SureTemp® Plus Electronic Thermometer**  
Provides accurate oral readings in about four to six seconds, 10- to 13-second pediatric axillary and rectal readings and 12- to 15-second adult axillary readings.

## The Welch Allyn Green Series™ 777 Integrated Wall System

### Recommended Systems:



**77791-MX**  
Includes GS 777 Wall Transformer (#77710), Coaxial Ophthalmoscope (#11720), MacroView Otoscope (#23810), Wall Aneroid (#7670-01), KleenSpec Dispenser (#52400-PF), and SureTemp Plus Thermometer (#01690-300).



**77791-2MPX**  
Includes GS 777 Wall Transformer (#77710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), Wall Aneroid (#7670-01), KleenSpec Dispenser (#52400-PF), and SureTemp Plus Thermometer (#01690-300).



**77791-2MPXPROBP**  
Includes GS 777 Wall Transformer (#77710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), KleenSpec Dispenser (#52400-PF), Connex ProBP Digital Blood Pressure Device (#34XFWT-B) and SureTemp Plus Thermometer (#01690-300).



**77794-2MPLXI**  
Includes GS 777 Wall Transformer (#77710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), KleenSpec Dispenser (#52400-PF), and Wall Mount for Spot Vital Signs LXi (#4500-62). Spot Vital Signs/Spot Vital Signs LXi ordered separately.

	WALL BOARD	GS 777 WALL TRANSFORMER	COAXIAL OPHTHALMOSCOPES	PANOPTIC™ OPHTHALMOSCOPES	DIAGNOSTIC MACROVIEW™ OTOSCOPES	PNEUMATIC OTOSCOPE	KLEENSPEC® DISPENSERS	WALL ANEROID	CONNEX® PROBP™ DEVICE	SURETEMP® PLUS THERMOMETER			
Part #		77710	11720	11810	11820	23810	23820	20200	52400-PF	52100-PF	7670-01	34XFWT-B	01690-300
<b>GS 777 Integrated Wall Systems</b>													
77798	•	•	•					•		•	•		
77791-1MPX	•	•		•		•			•		•		•
77791-2MPX	•	•			•	•			•		•		•
77791-MX	•	•	•			•			•		•		•
77792-1MP	•	•		•		•			•		•		
77792-2MP	•	•			•	•			•		•		
77792-M	•	•	•			•			•		•		
77795-2MPX	•	•			•		•		•		•		•
77795-MX	•	•	•				•		•		•		•
77796-M	•	•	•				•		•		•		
77797-2PX	•	•			•			•		•	•		•
77797-X	•	•	•					•		•	•		•
77791-MXNOBP	•	•	•			•			•				•
77792-MNOBP	•	•	•			•			•				
<b>GS 777 Integrated Wall Systems with Connex® ProBP™ Digital Blood Pressure Device (includes wall mounting bracket and device)</b>													
77791-MXPROBP	•	•	•			•			•			•	•
77791-2MPXPROBP	•	•			•	•			•			•	•
<b>GS 777 Integrated Wall Systems for Spot Vital Signs® or Spot Vital Signs® LXi Devices (includes device wall mounting bracket only, device(s) must be ordered separately)</b>													
77794-MSPOT	•	•	•			•			•				
77794-MLXI	•	•	•			•			•				
77794-2MPSPOT	•	•			•	•			•				
77794-2MPLXI	•	•			•	•			•				
<b>GS 777 Wall Transformer Sets</b>													
77510		•							•				
77910		•							•		•		
77710-71M		•	•			•							
77710-81M		•		•		•							
77710-82M		•			•	•							
77710-TAKX		•							•		•		•